

HOW TO USE PROBIOTIC BALL CLEANER

Machine (ball cleaning machine!) wash instruction:

- Make sure ball sets are NEVER mixed, nearly each set is unique
- Preparation (take all you need)
- Switch on/start the machine
- Once on full speed, spray gently 50cms above the balls (max. 2x triggering)
- Let the machine finish its cycle and wait 1 min.
- Wipe down the balls you took out with your hands with the BALL CLEANER MICROFIBER 4BILLIARD cloth or use gloves!

Manual wash instruction:

- Make sure ball sets are NEVER mixed, nearly each set is unique
- Preparation (take all you need)
- Put the ball set in a clean! bowl or on a platter
- Spray gently 50cms above the balls, they need to be well moist!
- Dry them individually with a clean FACILITY MICROFIBER 4BILLIARD cloth.
- Polish them with the BALL CLEANER MICROFIBER 4BILLIARD cloth or use gloves before bringing into play!

MAKE SURE CLEAN BALLS and CLEAN HANDS are used in play to ensure your total PROTECTION SHIELD our products bring 4BILLIARD!

4BILLIARD®

